

DATA SHEET

VETMAT KENNEL ASSEMBLY

Strong and long-lasting kennels, to face, without spoiling, the harsh conditions of veterinary hospitalization.

Sliding closure with self-closing, high security to avoid incorrect closings or leaks of animals.

Solid and seamless floor and door grid to avoid possible accidents.

The inner base of each kennel is tilted in a way that prevents liquid backflow and makes cleaning easier.

The anterior edges of the inner bases are elevated to prevent fluid leakage and cross infection.

The lower ckennel is double and has a removable divider, to convert it into two individual spacess, if required.

Base with wheels of high mobility and resistance, with brakes.

Modular system. Different combinations of kennels can be made, depending on the space and the specific needs of each hospital.



100% AISI stainless steel, anti-corrosion, anti-acid oxidation.

Steel plate thickness 1.2mm.

Door bars, combination of 8mm and 6mm bars.

Grid with combined 10mm and 4mm bars.

The wheels on the base are medical grade, high intensity use, with brakes.

Size:

Complete assembly: 1220*700*1570mm Inner grid: 690*545*10mm

Upper row: 2 x 610*610*700mm Doors upper kennels: 530*480*8mm

Lower kennel: 1220*820*700mm Doors lower kennel: 530*690*8mm

Collecting tray: 690*540*45mm

Item #	Description
07010027	VETMAT Kennel assembly





